

User manual

NAIL MACHINE START TO PRO



- | | |
|------------------------------------|-------------------|
| 1 Power switch (ON/OFF) | 5 Handpiece |
| 2 Forward / reverse switch | 6 Handpiece Stand |
| 3 Speed control knob for handpiece | 7 Control Unit |
| 4 Inlet connection for handpiece | 8 Indicator Lamp |

Prior to using the unit:

1. Read these instructions thoroughly and keep them safe for future reference.
2. Make sure the unit was not damaged in transit.
3. Store the carton and packaging material for any necessary future repair work.
4. Do not operate the unit:
 - if the power cord or plug are damaged
 - if the unit is damaged or has been dropped
 - if the unit does not function properly
5. Before using the unit, please make sure that the voltage indicated on the unit corresponds to main electricity (230V).
6. Install the unit on a flat stable surface. Set the unit in a clean dry area. Avoid placing the unit near any heating sources or steam cleaners.
7. Plug in the power cord of the machine to main plug.
8. Connect the handpiece to the inlet connection on the control unit.

Start up:

1. Turn on the on/off black switch on the machine. Machine should be at the lowest speed when starting. The indicator lamp on the control unit will be on.
2. Speed adjustment: select the desired speed by turning the variable speed control knob on the control unit.

Instrument changing:

1. To lock the tip: push the desired tip firmly into the hole at the rotating head of the handpiece until you hear two 'clicks'. This indicates a secured lock.
2. To remove the tip: turn off the device and gently pull out the tip by hand.

3. Only use burs with a 2.32mm to 2.35mm Ø shaft
Attention: Never use burs which are non-concentric, out of alignment or damaged. Only use burs with a maximum head diameter of 6mm. Never exceed the maximum speed recommended by the bur manufacturer.

Maintenance and Care:

1. Clean the controll unit with a soft moist cloth. Take care that no water enters the housing of the unit. After cleaning rub it dry using a soft dry cloth.
2. Keep your working surface free from dust to prevent any kind of dust infiltration into the handpiece.

Warning & Safety Instructions:

1. The unit is intended for use only by a qualified podiatrist and manicurist for podiatry and manicure treatments. It is not to be used for purposes other than podiatry and nail treatments. All warranties against damages to persons and property fail to exist in this case. Please contact Professionails NV for a detailed warranty policy.
2. Protect the device against water splashes. Should a liquid or a foreign object of any kind find its way into the unit, immediately switch off the unit and pull the mains plug. Consult your dealer for further information and technical assistance before resuming the operation of the unit.
3. Avoid electrical repairs.
4. Before performing any maintenance procedures always unplug the unit.
5. We recommend using only authorized repair services and replacing all parts with authentic spare parts.

Used symbols



Indoor use only



Confomité Européenne



Inform yourself about the local separate collection system for electronic products. Please act according to your local rules and do not dispose of your old products with your normal household waste.

Technical Data:

Input voltage:	220-240V - AC - 50Hz
Power input:	10W
Output voltage:	8-17V - DC
Power output:	6-9 VA
Maximum rotation speed:	25.000 RPM
Transport and storage temperature:	+ 5°C - +50°C
Operating temperature:	+10°C - +40°C
Insulation class:	

The consignment includes:

1 Control unit, 1 Hand piece, 1 Handpiece Stand and 3 Tips.